



Recombinant *Synechococcus elongatus* Orotidine 5'-phosphate decarboxylase (pyrF)

Product Code	CSB-MP656713FPY
Abbreviation	pyrF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31K20
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MTDLAARDRI IVPLDVPDVE AAITLIDRLP EAQFFKVGLE LRVATGSTVL TILKERQKKI FLDLKLHDIP NTMAGACRSA ARYGVDLLTI HATAGKPALE AAQAGAIAGA AAAGVQPPTL LAVTVLTSLS QQQLTEELHV PEAVPNYVQH LARQAEACGI GGCVCSPQEL VVLKTACSDD FVRVTPGVPR TWAAAGDQQR VMTPQQAIAA GATYLVVGRP ITAADDPAAA FQRVCEELTP
Source	Mammalian cell
Target Names	pyrF
Protein Names	Recommended name: Orotidine 5'-phosphate decarboxylase EC= 4.1.1.23 Alternative name(s): OMP decarboxylase Short name= OMPDCase Short name= OMPdecase
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.