



Recombinant *Synechococcus elongatus* D-fructose 1,6-bisphosphatase class 2/sedoheptulose 1,7-bisphosphatase (Synpcc7942_0505)

Product Code	CSB-MP656743FPY
Abbreviation	Synpcc7942_0505
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31QY2
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	MEKTIGLEII EVVEQAAIAS ARLMGKGEKN EADRVAVEAM RVRMNQVEML GRIVIGEGER DEAPMLYIGE EVGIYRDADK RAGVPAGKLV EIDIAVDPCE GTNLCAYGQP GSMVLAISE KGGLFAAPDF YMKKLAAPPA AKGKVDINKS ATENLKILSE CLDRAIDELV VVVMDRPRHK ELIQEIRQAG ARVRLISDGD VSAAISCGFA GTNTHALMGI GAAPEGVISA AAMRCLGGHF QGQLIYDPEV VKTGLIGESR ESNIARLQEM GITDPDRVYD ANELASGQEV LFAACGITPG LLMGVRRFFK GGARTQSLVI SSQSRTARFV DTVHMFDDVK TVSLR
Source	Mammalian cell
Target Names	Synpcc7942_0505
Protein Names	Recommended name: D-fructose 1,6-bisphosphatase class 2/sedoheptulose 1,7-bisphosphatase Short name= FB Pase class 2/SBPase EC= 3.1.3.11 EC= 3.1.3.37 Alternative name(s): D-fructose-1,6-bisphosphate 1-phosphohydrolase class 2
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.