



Recombinant *Synechococcus elongatus* Probable rRNA maturation factor (Synpcc7942_0402)

Product Code	CSB-BP656744FPY
Abbreviation	Synpcc7942_0402
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31R85
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	MALLDLTVQM ACSDRLPDLP WQDWLDRWLT DIQPQLPEDW QAPTYEACLR FVDDAEIQQL NRDYRQLDKP TDVLAFAALE DELALGFDPE EPLYLGDVII SVPTAQRQAQ GHALETENVW LSAHGLLHLL GWDHPNEEQ AAMLTQQRW LQLVAVPVP RWLDEVSQPP SSEP
Source	Baculovirus
Target Names	ybeY
Protein Names	Recommended name: Probable rRNA maturation factor
Expression Region	1-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.