



Recombinant Danio rerio Ashwin (im:6907446, si:dkeyp-118h3.6)

Product Code	CSB-YP656975DIL
Abbreviation	im:6907446, si:dkeyp-118h3.6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q32LR5
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MASHRTDRTK NPTSNQDVSK VDLLLLHPELL SQEFIQLMLQ ERNIAVSDPE DRDRLTGLYL QHVIPLPQRE LPRSRWGKRM EKSRPRLSSS STHSSSTDSG RKRPLIVFDG SSTSTGCVKL KKPSSSAPV TTDRLKPPVS SVSLTNPIRK LSGASTNCSS SNFSNRTPVS SSGAIPKSPS NHSNSSVHSN NATSKLKRTS PSKGEPDTAK DIKSPETKKK IQHITWP
Source	Yeast
Target Names	im:6907446
Protein Names	Recommended name: Ashwin
Expression Region	1-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.