



# Recombinant *Oryza sativa* subsp. japonica Probable chromo domain-containing protein LHP1 (LHP1)

<b>Product Code</b>	CSB-BP657039OFG
<b>Abbreviation</b>	LHP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q339W7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MARGKNHPDG EEEEAPAAAG EEEAPVEMDE EGEMEEEEEE EQGEGEGEER DEGESEEEWE DAEVEEGGEE SEQAAEEED SPLVVVEAEA AAPVVDGSP KLAEGYYEIE DIRRRRLRKG KLQYLKWRG WPESANTWEP LENLSACSDI IDAFEMRLQS PRPGRKRKRK ITTTPVAGSN PSHGKRGRPR LDKASHTRAP APEPKQLPCR TSCRATNCS SKTVAGLDAS GSVVRNQLAQ NIVQEGSSSV ISRTPCQELP LSIRLTDQQN EHHLVNGSSN SENLVKVPPS QGGQVTGAKK RKSGNVRRFE QNKPTQGQGE CGALVVAEDV GSTEGETGDK KKTEGCPNRV HITKIIKPVR FAAVNNDVQ QVSITFKALR SDGQEVMDDD KELKANNPLL LISYYEQCLR YNPTS
<b>Source</b>	Baculovirus
<b>Target Names</b>	LHP1
<b>Protein Names</b>	Recommended name: Probable chromo domain-containing protein LHP1 Alternative name(s): Protein LIKE HETEROCHROMATIN PROTEIN 1
<b>Expression Region</b>	1-415
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.