



Recombinant *Oryza sativa* subsp. japonica Probable chromo domain-containing protein LHP1 (LHP1)

Product Code	CSB-EP657039OFG
Abbreviation	LHP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q339W7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MARGKNHPDG EEEEAPAAAG EEEAPVEMDE EGEMEEEEEE EQGEGEGEER DEGEEEEWE DAEVEGGEE SEQAAEEED SPLVVVEAEA AAPVVDGSP KLAEGYYEIE DIRRRRLRKG KLQYLVKWRG WPESANTWEP LENLSACSDI IDAFEMRLQS PRPGRKRKRK ITTTPVAGSN PSHGKRGRPR LDKASHTRAP APEPKQLPCR TSCRATNCS SKTVAGLDAS GSVVRNQLAQ NIVQEGSSSV ISRTPCQELP LSIRLTDQQN EHHLVNGSSN SENLVKVPPS QGGQVTGAKK RKSGNVRRFE QNKPTQGQGE CGALVVAEDV GSTEGETGDK KKTEGCPNRV HITKIIKPVR FAAVNNDVQ QVSITFKALR SDGQEVMDDD KELKANNPLL LISYYEQCLR YNPTS
Source	E.coli
Target Names	LHP1
Protein Names	Recommended name: Probable chromo domain-containing protein LHP1 Alternative name(s): Protein LIKE HETEROCHROMATIN PROTEIN 1
Expression Region	1-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.