



Recombinant Arabidopsis thaliana Homeobox protein SHOOT MERISTEMLESS (STM)

Product Code	CSB-MP657125DOA
Abbreviation	STM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q38874
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MESGSNSTSC PMAFAGDNSD GPMCPMMMMM PPIMTSHQHH GHDHQHQQQE HDGYAYQSHH QQSSSLFLQS LAPPQGTKNK VASSSSPSSC APAYSLMEIH HNEIVAGGIN PCSSSSSSAS VKAKIMAHPH YHRLLAAYVN CQKVGAPPEV VARLEEACSS AAAAAASMGP TGCLGEDPGL DQFMEAYCEM LVKYEQELSK PFKEAMVFLQ RVECQFKSLS LSSPSSFSGY GETAIDRNNN GSSEEEVDMN NEFVDPQAED RELKGQLLRK YSGYLGSLKQ EFMKKRKKGK LPKEARQQLL DWWSRHYKWP YPSEQQKLAL AESTGLDQKQ INNWFINQRK RHWKPSEDMQ FVVM DATHPH HYFMDNVLGN PFPMDHISST ML
Source	Mammalian cell
Target Names	STM
Protein Names	Recommended name: Homeobox protein SHOOT MERISTEMLESS
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.