



# Recombinant Arabidopsis thaliana WD repeat-containing protein LWD2 (LWD2)

<b>Product Code</b>	CSB-MP657131DOA
<b>Abbreviation</b>	LWD2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q38960
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVTSSDQIQN GSEEQSKRSE IYTYEAPWQI YAMNWSIRRD KKYRLAITSL IEQYPNRVEI VQLDESNGEI RSDPNLCFEH PYPPTKTSFI PDKECQRPD LATSSDFLRL WRISDDESRV ELKSCSSDK NSEFSGPITS FDWNEAEP IGTSSIDTTC TIWDIEREVV DTQLIAHDKE VYDIAWGGVG VFASVSEDGS VRVFDLRDKE HSTIIYESGE PSTPLVRLSW NKQDPRYMAT VIMGSAKIVV LDIRFPALPV VELQRHQASV NAIAWAPHSS SHICSAGDDS QALIWDISSM GQHVEGGLDP ILAYTAGAEV EQLQWSSSQP DWVAIAFSNK LQILRV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	LWD2
<b>Protein Names</b>	Recommended name: WD repeat-containing protein LWD2 Alternative name(s): Protein ANTHOCYANIN 11-B Short name= AtAN11-B Protein LIGHT-REGULATED WD2 WD repeat-containing protein AN11-B
<b>Expression Region</b>	1-346
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.