



Recombinant *Synechococcus elongatus* 3-deoxy-manno-octulosonate cytidyltransferase (kdsB)

Product Code	CSB-MP657405FPY
Abbreviation	kdsB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31KU9
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MVRILAVIPA RYASERLPGK VLLPIAGRPM LQWVYEATIA SNVFDQVAIA TEDPRVVEAA AAFGAEAILT SADLASGTDR VAEASLHFPD CKVIANVQGD QPFVTPGLLQ ALVSPYRAGE LPEMTTVGGP YDPAQDADDP NTVKVVCDQR GNALYFSRSA IPYPRTVVHD LPVYHHFGLY AFRRDFLAQY RQLPPTPLER CESLEQLRVL EQGYRIRVVP CADKVIEVNT ADDLERANAW ASQR
Source	Mammalian cell
Target Names	kdsB
Protein Names	Recommended name: 3-deoxy-manno-octulosonate cytidyltransferase EC=2.7.7.38 Alternative name(s): CMP-2-keto-3-deoxyoctulosonic acid synthase Short name= CKS Short name= CMP-KDO synthase
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.