



Recombinant *Synechococcus elongatus* Phosphomethylpyrimidine synthase (thiC)

Product Code	CSB-BP657423FPY
Abbreviation	thiC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31P93
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	MRSDWIAPRR GQANVTQMHY ARQGVITEEM DFVARRENLP ADLIRDEVAR GRMIIPANIN HTNLEPMAIG IASKCKVNAN IGASPNASNI DEEVEKLLA VKYGADTVMD LSTGGGNLDE IRTAINASP VPIGTVPVYQ ALESVHGRIE KLSADDFLHV IEKHCEQGVD YQTIHAGLLI EHLPKVKSRI TGIVSRGGGI IAQWMLYHHK QNPLYTHFRD IIEIFKRYDC SFSLGDSL RP GCLHDASDDA QLSELKTLGQ LTRVAWEHDV QVMVEGPGHV PMDQIEFNVR KQMEECSEAP FYVLGPLVTD IAPGYDHITS AIGAAMAGWY GTAMLCYVTP KEHLGLPNAE DVRNGLIAYK IAAHAADIAR HRPGARDRDD ELSRARYAFD WNKQFDLSLD PERAREYHDE TLPADIYKTA EFCSMCGPKH CPMQTKITEE DLTELEKFLE KDSALA
Source	Baculovirus
Target Names	thiC
Protein Names	Recommended name: Phosphomethylpyrimidine synthase EC= 4.1.99.17 Alternative name(s): Hydroxymethylpyrimidine phosphate synthase Short name= HMP-P synthase Short name= HMP-phosphate synthase Short name= HMPP synthase
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.