



Recombinant *Synechococcus elongatus* 4-hydroxy-3-methylbut-2-enyl diphosphate reductase (ispH)

Product Code	CSB-EP657436FPY-B
Abbreviation	ispH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31S64
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MDTKAFKRAL HQSDRYNRKG FGKTTDVSGA LESAYQSDLI QSLRQNGYRL QRGEITIRLA EAFGFCWGVE RAVAIAYETR QHFPQERIWI TNEIIHNPSV NQHLREMSVE FIPCERGEKD FSVVDRGDVV ILPAFGASVQ EMQLLDEKGC HIVDTTCPWV SKVWNTVEKH KRGHTSIIH GKYNHEETVA TSSFAETYLV VLNLEQAQYV CDYILNGGDR DEFMTRFGKA CSAGFDPDRD LERIGIANQT TMLKSETEAI GKLFFERTLLK KYGPQALNDH FLAFNTICDA TQERQDAMFQ LVEEPLDLIV VIGGFNSSNT THLQEIAIER QIPSFHIDAA ERIGPGNRIE HKPLHTDLTT TEPWLPAGPL TIGITSGAST PDKVVEDVIE RLFDLQRS
Source	E.coli
Target Names	ispH
Protein Names	Recommended name: 4-hydroxy-3-methylbut-2-enyl diphosphate reductase EC= 1.17.1.2
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.