



Recombinant *Oryza sativa* subsp. *japonica* Uncharacterized protein ycf73 (ycf73-A)

Product Code	CSB-EP657786OFG-B
Abbreviation	ycf73-A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q37066
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MTKDETLVVF TLVSSVSVF LFGILLFMVL ISATRDFRER TSKSLVKIMI WAGIVVITFA IAVRIYPIFI FLLKERIKPL VEALYDKLPW IWEVLSRYW DRLIDFLDRY LWACAQRIQT GIRKQKGEFV VTFSCRVKKR LYARAIEVGI HLSLLSNLFW ILKTTLAVGY RLLWVLYYII SFEGFLGSFR LYLVEYFGFYC LLFSGKWLRT SEDRGERQAQ ISGILLRGML IECAFSVLCL EEDSNLHAL
Source	E.coli
Target Names	ycf73-A
Protein Names	Recommended name: Uncharacterized protein ycf73 Alternative name(s): ORF249
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.