



Recombinant Arabidopsis thaliana Probable xyloglucan endotransglucosylase/hydrolase protein 15 (XTH15)

Product Code	CSB-EP657822DOA
Abbreviation	XTH15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q38911
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	SNFFD EFDLTWGDHR GKIFNGGNML SLSLDQVSGS GFKSKKEYLF GRIDMQLKLV AGNSAGTVTA YYLSSQGATH DEIDFEFLGN ETGKPYVLHT NVFAQGKGDREQQFYLWFDPTKNFHTYSIVWRPQHIIFLV DNLPIRVFNN AEKLGVPFPK SQPMRIYSSL WNADDWATRG GLVKTDWSKA PFTAYYRGFN AACTASSGC DPKFKSSFGEKLVATELNAYGRRRLRWV QKYFMIYNYC SDLKRFPRGF PPECKKSRV
Source	E.coli
Target Names	XTH15
Protein Names	Recommended name: Probable xyloglucan endotransglucosylase/hydrolase protein 15 Short name= At-XTH15 Short name= XTH-15 EC= 2.4.1.207
Expression Region	26-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.