



Recombinant Xenopus laevis Tumor necrosis factor receptor type 1-associated DEATH domain protein (tradd)

Product Code	CSB-EP658365XBE-B
Abbreviation	tradd
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q32NG6
Product Type	Recombinant Protein
Immunogen Species	Xenopus laevis (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MAAGPSMWVG SVYLYIKSDT VPLPGKYTHQ KALIYEALRS AISESTRGCR DSIEILKIHS SDQQLILYLK FCGLEPCQRF LKDYKECKVQ MGIQNKLKNC FSVEGLPIFT ELKIDTGEID SLLEKEEQCL KYISQMKPTI QKDELAEID ERLKSILKDS PSALDSELSL QNSCQCSSLPL SLHSNRSYHI EGSTFHFQGE EFVDRPLTSA HIQHFAKSVG KNWKPVGRSL GKTCRALNDT AIENLAYEFD RDGRYEQAYQ LLRLFKDSEG KKATVQRLVQ ALEENGLNSI ALDLLSLNEN GLK
Source	E.coli
Target Names	tradd
Protein Names	Recommended name: Tumor necrosis factor receptor type 1-associated DEATH domain protein Short name= TNFR1-associated DEATH domain protein Alternative name(s): TNFRSF1A-associated via death domain
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.