



Recombinant Xenopus laevis Tumor necrosis factor receptor type 1-associated DEATH domain protein (tradd)

| | |
|--------------------------|---|
| Product Code | CSB-EP658365XBE |
| Abbreviation | tradd |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q32NG6 |
| Product Type | Recombinant Protein |
| Immunogen Species | Xenopus laevis (African clawed frog) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MAAGPSMWVG SVYLYIKSDT VPLPGKYTHQ KALIYEALRS AISESTRGCR DSIEILKIHS SDQQLILYLK FCGLEPCQRF LKDYKECKVQ MGIQNKLKNC FSVEGLPIFT ELKIDTGEID SLLEKEEQCL KYISQMKPTI QKDELAEID ERLKSILKDS PSALDSELSL QNSCQCSSLPL SLHSNRSYHI EGSTFHFQGE EFVDRPLTSA HIQHFAKSVG KNWKPVGRSL GKTCRALNDT AIENLAYEFD RDGRYEQAYQ LLRLFKDSEG KKATVQRLVQ ALEENGLNSI ALDLLSLNEN GLK |
| Source | E.coli |
| Target Names | tradd |
| Protein Names | Recommended name: Tumor necrosis factor receptor type 1-associated DEATH domain protein Short name= TNFR1-associated DEATH domain protein Alternative name(s): TNFRSF1A-associated via death domain |
| Expression Region | 1-303 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.