



Recombinant Arabidopsis thaliana Agamous-like MADS-box protein AGL15 (AGL15)

Product Code	CSB-EP658511DOA
Abbreviation	AGL15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q38847
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MGRGKIEIKR IENANSRQVT FSKRRSGLLK KARELSVLCD AEVAVIVFSK SGKLFEYSST GMKQTLSTRYG NHQSSSASKA EEDCAEVDIL KDQLSKLQEK HLQLQGKGLN PLTFKELQSL EQQLYHALIT VRERKERLLT NQLEESRLKE QRAELENETL RRQVQELRSF LPSFTHYVPS YIKCFAIDPK NALINHDSKC SLQNTDSDTT LQLGLPGEAH DRRTNEGERE SPSSDSVTTN TSSETAERGD QSSLANSPPE AKRQRFSV
Source	E.coli
Target Names	AGL15
Protein Names	Recommended name: Agamous-like MADS-box protein AGL15
Expression Region	1-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.