



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At3g26630, chloroplastic (PCMP-A6)

Product Code	CSB-BP658521DOA
Abbreviation	PCMP-A6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q38959
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	S PEASYFLRTC SNFSQLKQIH TKIHKHNLTN DQLLVRQLIS VSSSFGETQY ASLVFNQLQS PSTFTWNLMI RLSLVNHKPR EALLFILMM ISHQSQFDKF TFPFVIKACL ASSSIRLGTQ VHGLAIKAGF FNDVFFQNTL MDLYFKCGKP DSGRKVFDMK PGRSIVSWTT MLYGLVSNSQ LDSAEIVFNQ MPMRNVVSWT AMITAYVKNR RPDEAFQLFR RMQVDDVKNP EFTIVNLLQA STQLGSLSMG RWVHDYAHKN GFVLDCFLGT ALIDMYSKCG SLQDARKVFD VMQGKSLATW NSMITS LGVH GCGEEALSLF EEMEEEASVE PDAITVGVVL SACANTGNVK DGLRYFTRMI QVYGISPIRE HNACMIQLLE QALEVEKASN LVESMDSDPD FNSSFGNEYT DGMNETNETP SQHQIMFTKW DTGRF
Source	Baculovirus
Target Names	PCMP-A6
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At3g26630, chloroplastic
Expression Region	20-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.