



Recombinant Human DNA-directed RNA polymerase I subunit RPA43 (TWISTNB)

Product Code	CSB-YP658940HU
Abbreviation	TWISTNB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3B726
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAGCSEAPR PAAASDGSLV GQAGVLPCL LPTYAAACAL VNSRYSCLVA GPHQRHIALS PRYLNRKRTG IREQLDAELL RYSESLLGVP IAYDNIKVVG ELGDIYDDQG HIHLNIEADF VIFCPEPGQK LMGIVNKVSS SHIGCLVHGC FNASIPKPEQ LSAEQWQTME INMGDELEFE VFRLDSDAAG VFCIRGKLNI TSLQFKRSEV SEEVTEGTE EAAKKPKKKK KKKDPETYEYV DSGTTKLADD ADDTPMEEES LQNTNNANGI WEEEPKKKKK KKKHQEVQDQ DPVFQGS DSS GYQSDHKKKK KKRKHSEEA FTPPLKCS PK RKGKSNFL
Source	Yeast
Target Names	TWISTNB
Protein Names	Recommended name: DNA-directed RNA polymerase I subunit RPA43 Alternative name(s): Twist neighbor protein
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.