



# Recombinant Rat Glycogen synthase kinase-3 beta (Gsk3b)

<b>Product Code</b>	CSB-EP009963RA-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P18266
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSGRPRTTSF AESCKPVQQP SAFGSMKVS R DKDGSKVTTV VATPGQGPDR PQEVSYTDTK VINGNSFGVV YQAKLCDSGE LVAIKKVLQD KRFRNRELQI MRKLDHCNIV RLRYFFYSSG EKKDEVYLN L VLDYVPETVY RVARHYSRAK QTLPLVIYVKL YMYQLFRSLA YIHSFGICHR DIKPQNLLLD PDTAVLKLCD FGS AKQLVRG EPNVSYICSR YYRAPELIFG ATDYTSSIDM WSAGCVLAEL LLGQPIFPGD SGVDQLVEII KVLGTP TREQ IREMNP NYTE FKFPQIKAHP WTKVFRP RTP PEAIALCSRL LEYTP TARLT PLEACAHSFF DEL RDPNVKL PNGRDT PALF NFFTQELSSN PPLATILIPP HARIQAAASP PANATAASDT NAGDRGQTNN AASASASNST
<b>Source</b>	E.coli
<b>Target Names</b>	Gsk3b
<b>Protein Names</b>	Recommended name: Glycogen synthase kinase-3 beta Short name= GSK-3 beta EC= 2.7.11.26 Alternative name(s): Factor A Short name= FA Serine/threonine-protein kinase GSK3B EC= 2.7.11.1
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.