



Recombinant Human ADP-ribosylation factor-like protein 13B (ARL13B)

Product Code	CSB-MP659639HU
Abbreviation	ARL13B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3SXY8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFSLMASCCG WFKRWREPVR KVTLLMVGLD NAGKTATAKG IQGEYPEDVA PTVGFSKINL RQGKFEVTIF DLGGGIRIRG IWKNYYAESY GVIFVVDSSD EERMEETKEA MSEMLRHPRI SGKPIVLAN KQDKEGALGE ADVIECLSLE KLVNEHKCLC QIEPCSAISG YGKKIDKSIK KGLYWLLHVI ARDFDALNER IQKETTEQRA LEEQEKQERA ERVRKLREER KQNEQEQAEL DGTSGLAELD PEPTNPFQPI ASVIIENEGK LEREKKNQKM EKDSDGCHLK HKMEHEQIET QGQVNHNGQK NNEFGLVENY KEALTQQLKN EDETDRPSLE SANGKKKTKK LRMKNRHRVE PLNIDDCAPE SPTPPPPPPP VGWGTPKVTR LPKLEPLGET HHNDFYRKPL PPLAVPQRPN SDAHDDVIS
Source	Mammalian cell
Target Names	ARL13B
Protein Names	Recommended name: ADP-ribosylation factor-like protein 13B Alternative name(s): ADP-ribosylation factor-like protein 2-like 1 Short name= ARL2-like protein 1
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.