



Recombinant Mouse Eukaryotic translation initiation factor 4E type 1B (Eif4e1b)

Product Code	CSB-EP659749MO
Abbreviation	Eif4e1b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3UTA9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNKVEGGGHK EEVVVKEKEV VKEKPSEATA EGVQAGEAKD LPGSLKTQRR KAHREHPPEV LSKLHPLQYR WVLWFFKNDR SRAWQDNLQL VTKFNTVEDF WAVYSHIKLA SKLSSGCDYA LFKEGILPMW EDNRNKQGGR WLLSIDKQLR HFELDRLWLE TLLCLVGNCF EEYSREVCGA VVNIRTKRDK IALWTSEAED KAGVMQIGQI YKERLGISTK TIIGYQAHAD TAAKSNLAN KFVV
Source	E.coli
Target Names	Eif4e1b
Protein Names	Recommended name: Eukaryotic translation initiation factor 4E type 1B Alternative name(s): Oocyte-specific eukaryotic translation initiation factor 4E-like
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.