



Recombinant Solanum tuberosum Patatin-05 (pat1-k1)

Product Code	CSB-MP659800FIG
Abbreviation	pat1-k1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3YJT5
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	TLGEMVT VLSIDGGGIK GIIPGIILEF LEGQLQKMDN NADARLADYF DVIGGTSTGG LLTAMITTPN ENNRPFAAAK DIVPFYFQHG PHIFNSSTGQ FFGPKYDGKY LMQVLQEKLK ETRVHQALTE VAISSFDIKT NKPVIFTKSN LAKSPELDAK MSDICYSTAA APTYFPPHYF ATNTSNGDKY EFNLVDGAVA TVADPALLSV SVATRRAEED PAFASIRSLN YKQLLLLSLG TGTNSEFDKT HTAQETAKWG ALQWMLVIQQ MTEAASSYMT DYYLSTVFDQ LHSQNNYLRV QENALTGTTT KADDASEANM ELLVQVGENL LKKPVSKDNP ETYEEALKRF AKLLSDRKKF RANKASY
Source	Mammalian cell
Target Names	pat1-k1
Protein Names	Recommended name: Patatin-05 EC= 3.1.1.- Alternative name(s): Patatin-1-Kuras 1
Expression Region	24-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.