



Recombinant Arabidopsis thaliana Probable FBD-associated F-box protein At5g38565 (At5g38565)

Product Code	CSB-MP660356DOA
Abbreviation	At5g38565
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3E8L3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDIFNGLPDD VLVKILSFVP TKVAVSTSIL SKRWEFLWMW LPRLDGSPK TDLLAFLNTC QYDKEEEVGL VDFIDKKLPL HRAPFLAMYV CWKFLMNVGM KTYFHVRHLE TRQREGKSLQ DILSICPVL DLSVICSVHQ DVKEFTIIVP SLQSLTLFIE NCEVFDGYVI DTPLKYLKLE DVHEEEHYCL LKKMPKLREA YVDVQLDDLK SLIGSITSVK RLNHMFREKS KILGQLLKDS PNLRLVNIK VQGHVTLSTG VDCWNQPISV PECLLESQI FNLSHYFGKQ QDLDFVYIIL KNACHLKTAT ILADEPEHLV PNLKELTLSP RASSTCQLSI RCGLGERS
Source	Mammalian cell
Target Names	At5g38565
Protein Names	Recommended name: Probable FBD-associated F-box protein At5g38565
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.