



Recombinant Arabidopsis thaliana F-box protein At1g67130 (At1g67130)

Product Code	CSB-MP660366DOA
Abbreviation	At1g67130
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3ECH0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MNRGETLDSI PTDLILDILS RLPTKSIARF HCVSKLWSSM LASQDFTRLF VNRSSSNPRL LLGIVRGGEW SFYSSPQPKN PYEKSSLEVA ADFHMKLSEI KQHQHHGTSY AFGLIYLRTV WMSKEGNYEV CVICKPSTGQ YYAILPPQPS GIFGDLGFDP IGKQFKVLVL NILDNHYILT LGTENDKWRS IQSSLRYRPC GQSPICINGV LYYIAYDTQD SSNDVIGCFD VRFEKFKFFH VNHDMVKCFF ELINYKDKLG GFPLELSMWV LEDLEKEEWS KYAYTLKPDN NVVKVNYNLS FVGVTTRGEI VLAKMYACKP FYVFYYPEN NTLLSVEIQG VGEDSEWFSN YQRVYVFVEH VKDLHQFDNT KTSINLPEQK RKPTSISISS KYDDHVRTTL ISSRKNQQVT TVSRPQQDRR TTNKFSALCL LVMKNSQVSK PDCFSWSSAT MLYL
Source	Mammalian cell
Target Names	At1g67130
Protein Names	Recommended name: F-box protein At1g67130
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.