



Recombinant Human Zinc finger protein SNAI3 (SNAI3)

Product Code	CSB-BP660752HU
Abbreviation	SNAI3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3KNW1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPRSFLVKTH SSHRVPNYRR LETQREINGA CSACGGLVVP LLPRDKEAPS VPGDLPQPWD RSSAVACISL PLLPRIIEAL GASGLDALEV SEVDPRASRA AIVPLKDSLNL HLNLPPLLVL PTRWSPTLGP DRHGAPEKLL GAERMPRAPG GFECFHCHKP YHTLAGLARH RQLHCHLQVG RVFTCKYCDK EYTSLGALKM HIRTHLTPCT CKICGKAFSR PWLLQGHVRT HTGEKPYACS HCSRAFADRS NLRAHLQTHS DAKKYRCRRC TKTFSRMSLL ARHEESGCCP GP
Source	Baculovirus
Target Names	SNAI3
Protein Names	Recommended name: Zinc finger protein SNAI3 Alternative name(s): Protein snail homolog 3 Zinc finger protein 293
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.