



Recombinant Human Ankyrin repeat domain-containing protein 55 (ANKRD55)

Product Code	CSB-EP660753HU-B
Abbreviation	ANKRD55
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3KP44
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMRQATMDFS TPSVFDQQRG DSSEEVDLTM VYQAASNGDV NALTAVIRED PSILECCDSE GCTPLMHAVS GRQADTVKLL LKMGANINMQ DAYGRTSLCL ATYLGWLEGC VLLLRNGAKH NIPDKNGRLP LHAATAEPDM RLLTVLLQQS NISEINHQDN EGMTPLHWAA FHNQPQHTQM LLKKGADPTL VDKDFKTALH WAVQSGNRIL CSIILSHHQG PSIINYDDES GKTCVHIAAA AGFSDIIHEL ARVPECNLQA LDVDDRTPLH WAAAAGKAEC VQSLLELGMD SNLRDINEST PLAYALYCGH TACVKLLSQE SRTEPTRPPP SQSSRPQKKE RRFNVLNQIF CKNKKEEQRA HQKDPSRDRY REEDTSEVND IITTFDSIVG TNCQEQPGDQ VAMVEFKKKT SDNSKYLLPE KKPLARKGLP PIRTQSLPPI TLGNNFLTAS HRATSHAGLS SAPHHMAQRS QKSRSEQDLL NNRTGCQMLL DNPWKSDSNQ VFSYKVVTVS SSDKLLDRLL SVRPGHQEVS VPPHLRHLHN PSSGQNFQHL SPNRHKIRD L PFTRNNLAPL PDQKFLSGEP LRTNRVLPAL PSQRRHSTAA EESEHSANPT SDEN
Source	E.coli
Target Names	ANKRD55
Protein Names	Recommended name: Ankyrin repeat domain-containing protein 55
Expression Region	1-614
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.