



Recombinant Human Testis-specific protein 10-interacting protein (TSGA10IP)

Product Code	CSB-MP660950HU
Abbreviation	TSGA10IP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3SY00
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGQDTDMLNT YQQLVRTPSV RPGQDVRLQA PGTRTGLLKL LSTVSQDKQG CLGSGDGVPN QDLQQRSQSS RQTAKKDRKP RGQSKKGQGS EESEDHFPLL PRKPSFPFQW AWESIATDVR AVLQPSSPTP GHQALPMPSS FSQRQSRKRS TANLPEAHGC CWKTEAQLNK ARQQLGAWGG VSIPTGKGEL GSEPPSGLQL PGRRPGSGSA SDKQVQLQSL GAEEAERGLS SGVLPQRPRR GSISEEEQFS EATEEAEEGE HRTPCRRRAG CQRKGQISGE EASDEGEVQG QSQGSSPSFN NLRRRQWRKT RAKELQGPWD LEKLHRQLQR DLDCGPQKLP WKTLRAAFQA SKRNGKAYAS GYDETFVSAN LPNRTFHKRQ EATRSLQAW ERQRQEERQQ AELRRARTQH VQRQVAHCLA AYAPRGSRGF GAAQRKLEEL RRQERQRF AE YQAE LQGIQH RVQARPFLFQ QAMQANARLT VTRRFSQVLS ALGLDEEQLL SEAGKVDREG TPRKPRSHRS VGV RMEHSPQ RPPRTEPTGS QPDRHYNPSL DPECSP
Source	Mammalian cell
Target Names	TSGA10IP
Protein Names	Recommended name: Testis-specific protein 10-interacting protein Alternative name(s): Tsga10-interacting protein
Expression Region	1-556
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.