



# Recombinant Human Testis-specific protein 10-interacting protein (TSGA10IP)

<b>Product Code</b>	CSB-YP660950HU
<b>Abbreviation</b>	TSGA10IP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3SY00
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGQDTDMLNT YQQLVRTPSV RPGQDVRLQA PGTRTGLLKL LSTVSQDKQG CLGSGDGVPN QDLQQRSQSS RQTAKKDRKP RGQSKKGQGS EESEDHFPLL PRKPSFPFQW AWESIATDVR AVLQPSSPTP GHQALPMPSS FSQRQSRRKS TANLPEAHGC CWKTEAQNLIK ARQQLGAWGG VSIPTGKGEL GSEPPSGLQL PGRRPGSGSA SDKQVQLQSL GAEEAERGLS SGVLPQRPRR GSISEEEQFS EATEEAEEGE HRTPCRRRAG CQRKGQISGE EASDEGEVQG QSQGSSPSFN NLRRRQWRKT RAKELQGPWD LEKLHRQLQR DLDCGPQKLP WKTLRAAFQA SKRNGKAYAS GYDETFVSAN LPNRTFHKRQ EATRSLQAW ERQRQEERQQ AELRRARTQH VQRQVAHCLA AYAPRGSRGP GAAQRKLEEL RRQERQRFAE YQAEIQGIQH RVQARPFQ QAMQANARLT VTRRFSQVLS ALGLDDEEQLL SEAGKVDREG TPRKPRSHRS VGVMEHSPQ RPPRTEPTGS QPDRHYNPSL DPECSP
<b>Source</b>	Yeast
<b>Target Names</b>	TSGA10IP
<b>Protein Names</b>	Recommended name: Testis-specific protein 10-interacting protein Alternative name(s): Tsga10-interacting protein
<b>Expression Region</b>	1-556
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.