



# Recombinant Bovine Complement factor D (CFD)

<b>Product Code</b>	CSB-BP660982BO
<b>Abbreviation</b>	CFD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3T0A3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ILRG QEAPSHSRPY MASVQVNGKH VCGGFLIAEQ WVMSAAHCLE DVADGKVQVL LGAHSLSQPE PSKRLYDVLV VVPHPGSRTE TIDHLLLLLQ LSEKAVLGPA VQLLPWQRED RDVAAGTLCD VAGWGVVSHT GRKPDRLQHL LLPVLDRTC NLRTYHDGTI TERMMCAESN RRDTCKGDSG GPLVCGSVAE GVTSGSRIC GNHKKPGIYT RLASYVAWID GVMAEGAAA
<b>Source</b>	Baculovirus
<b>Target Names</b>	CFD
<b>Protein Names</b>	Recommended name: Complement factor D EC= 3.4.21.46 Alternative name(s): Adipsin C3 convertase activator Properdin factor D
<b>Expression Region</b>	27-259
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a member of the trypsin family of peptidases. The encoded protein is a component of the alternative complement pathway best known for its role in humoral suppression of infectious agents. This protein is also a serine protease that is secreted by adipocytes into the bloodstream. Finally, the encoded protein has a high level of expression in fat, suggesting a role for adipose tissue in immune system biology.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.