



Recombinant Mouse Myosin regulatory light chain 12B (Myl12b)

Product Code	CSB-EP661006MO-B
Abbreviation	Myl12b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3THE2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSSKKAKTKT TKKRPQRATS NVFAMFDQSQ IQEFKEAFNM IDQNRDGFID KEDLHDMLAS LGKNPTDAYL DAMMNEAPGP INFTMFLTMF GEKLNQTDPE DVIRNAFACF DEEATGTIQE DYLLRELLTTM GDRFTDEEVD ELYREAPIDK KGNFNIEFT RILKHGAKDK DD
Source	E.coli
Target Names	Myl12b
Protein Names	Recommended name: Myosin regulatory light chain 12B Alternative name(s): Myosin regulatory light chain 2-B, smooth muscle isoform Myosin regulatory light chain 20 kDa Short name= MLC20 Myosin regulatory light chain MRLC2
Expression Region	1-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.