



Recombinant Mouse Tubulin alpha chain-like 3 (Tubal3)

Product Code	CSB-EP661063MO-B
Abbreviation	Tubal3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3UX10
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MRECLSIHIG QAGVQIGDAC WELYCLEHGI QPDGFILDHQ HDNLENPKVE HMNASLDFFF HETRAGKHVP RTLFMDLEPT VIDGIRVGRY HSLFHPEQLV NGKEDAANTY ARGYSVGSE VIELVLERIR KLAEQCSGLQ GFLIYRSFGG GTGSGFTSLL MERLSVEYCK KIKLEFSVYP SPRISTAVVE PYNAILTTHS TIEYSDCAFM VDNEALYDIC QHKLGIERPS YASINRLIAQ VSSSITASLR FEGPLNVDLI EFQTNLVPYP RIHPITALA PIISAEKAYQ EQLSVSDVTA SCFEVSNQLV KCDPRLGKYM ACCLLYRGDV VPKDVNEAIA AMKSRTSVQF VDWCPTGFKV GINYQPPAVV PGGDLARVQR AVCMLSNTTA IVEAWARLDH KFDLMYAKKA FLHWYITEGM ELGEFVEARE DLAALEKDYE EVGLSF
Source	E.coli
Target Names	Tubal3
Protein Names	Recommended name: Tubulin alpha chain-like 3
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.