



# Recombinant Mouse Pyrin domain-containing protein 3 (Pydc3)

<b>Product Code</b>	CSB-EP661077MO-B
<b>Abbreviation</b>	Pydc3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3V3Q4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTGEMVNYK QIVLLSGLEY MNDYNFRALK SLLNHDLKLT KNMQDDYDRI KIADLMEEKF PEDAGLSKLI EVCEDIPELA ARVDILRKEM EKVKNKTKIK SESSPLPLTS SLMEAWVEVP AMVTASSEES KDTIPESPD TTTQFLEEKP KFPLLSATST SQAEGEPLTP QRFPTTASSS LQTPLEPTEI LSTILATSQG SSAPYSTCDK SSRVPPVTVS SSLQTIQTCQ ATSTLPCSHH TSLESPKTEP SSVQATQMTQ AIKASGHNCP QVPASAVSSS FIKPQVTPAM LLSGVQTPLM PQATVPSRAQ TFQLTPAKMT SGCNSPQMSA ATVYSSYSNP QVTPVHSSVQ ILQMNLAAMT IGCNSPHVSA ATVSSPYNNP WVTPATFPRN AQTLLQLYPAA MAHACNSPQV SAATISSSYN NTPQVSSVTV PRSFPAMSL S PAKPLKAKIG SHLGATDQLV KDHQFEDNEM QNPQSGGLGTG LSDQPRLSRT QSRTEKRRRK QEEQPLATKM KRSLSPPTTARRVNPMP TARRVNPMP WGQLKLTRE AEGLVQRTGN KLSSETMFLA MLALIAMQSS SEICNYDKGG HSGTSSEK
<b>Source</b>	E.coli
<b>Target Names</b>	Pydc3
<b>Protein Names</b>	Recommended name: Pyrin domain-containing protein 3 Alternative name(s): Interferon-activable protein 204-like Short name= lfi-204-like
<b>Expression Region</b>	1-588
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.