



Recombinant Human Protein FAM161A (FAM161A)

Product Code	CSB-BP661553HU
Abbreviation	FAM161A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3B820
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MATSHRVAKL VASSLQTPVN PITGARVAQY EREDPLKALA AAEAILEDDEE EEKVAQPAGA SADLNTSFSG VDEHAPISYE DFVNFDPDIHH SNEEYFKKVE ELKAAHIETM AKLEKMYQDK LHLKEVQPVV IREDSLSDSS RSVSEKNSYH PVSLMTSFSE PDLGQSSSLY VSSSEEELPN LEKEYPRKNR MMTYAKELIN NMWTFDFCVED YIRCKDTGFH AAEKRRKKRK EWWPTITVPE PFQMMIREQK KKEESMKSXS DIEMVHKALK KQEEDPEYKK KFRANVPAS VFLPLYHDLV KQKEERRRSL KEKSKEALLA SQKPFKFIAR EEQKRAAREK QLRDFLKYKK KTNRFKARPI PRSTYGSTTN DKLKEEELYR NLRTQLRAQE HLQNSSPLPC RSACGCRNPR CPEQAVKLKC KHKVRCPTPD FEDLPERYQK HLSEHKSPKL LTVCKPFDLH ASPHASIKRE KILADIEADE ENLKETRWPY LSPRRKSPVR CAGVNPVPCN CNPPVPTVSS RGREQAVRKS EKERMREYQR ELEEREELK KRPLLFERVA QKNARMAAEK HYSNTLKALG ISDEFVSKKG QSGKVLEYFN NQETKSVTED KESFNEEEEKI EERENGEENY FIDTNSQDSY KEKDEANEE EEKSVEESH
Source	Baculovirus
Target Names	FAM161A
Protein Names	Recommended name: Protein FAM161A
Expression Region	1-660
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.