



# Recombinant Human Protein FAM161A (FAM161A)

<b>Product Code</b>	CSB-EP661553HU-B
<b>Abbreviation</b>	FAM161A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3B820
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MATSHRVAKL VASSLQTPVN PITGARVAQY EREDPLKALA AAEAILEDDEE EEKVAQPAGA SADLNTSFSG VDEHAPISYE DFVNFPDIHH SNEEYFKKVE ELKAAHIETM AKLEKMYQDK LHLKEVQPVV IREDSLSDSS RSVSEKNSYH PVSLMTSFSE PDLGQSSSLY VSSSEEELPN LEKEYPRKNR MMTYAKELIN NMWTFDFCVD YIRCKDTGFH AAEKRRKKRK EWWPTITVPE PFQMMIREQK KKEESMKS KS DIEMVHKALK KQEEDPEYKK KFRANVPAS VFLPLYHDLV KQKEERRRSL KEKSKEALLA SQKPFKFIAR EEQKRAAREK QLRDFLKYKK KTNRFKARPI PRSTYGSTTN DKLKEEELYR NLRTQLRAQE HLQNSSPLPC RSACGCRNPR CPEQAVKLKC KHKVRCPTPD FEDLPERYQK HLSEHKSPKL LTVCKPFDLH ASPHASIKRE KILADIEADE ENLKETRWPY LSPRRKSPVR CAGVNPVPCN CNPPVPTVSS RGREQAVRKS EKERMREYQR ELEEREELK KRPLLFERVA QKNARMAAEK HYSNTLKALG ISDEFVSKKG QSGKVLEYFN NQETKSVTED KESFNEEEEKI EERENGEENY FIDTNSQDSY KEKDEANES EEKSVEESH
<b>Source</b>	E.coli
<b>Target Names</b>	FAM161A
<b>Protein Names</b>	Recommended name: Protein FAM161A
<b>Expression Region</b>	1-660
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.