



Recombinant *Phalaenopsis aphrodite* subsp. *formosana* Maturase K (matK)

Product Code	CSB-MP661565PDAO
Abbreviation	matK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3BAR2
Product Type	Recombinant Protein
Immunogen Species	<i>Phalaenopsis aphrodite</i> subsp. <i>formosana</i> (Moth orchid)
Purity	>85% (SDS-PAGE)
Sequence	MYQQKYLISS VSDSNQNGFS LHKNSFSYHF YSQMVSEGFG VILEISFSSR LVSSLEEKRI PKSQNLRSIH SIFPFLEDKL SHLNYVSDLL IPYPIHMEIL VQILQCWIKD VPSLHLLRFI FHEYHNLNSL ITSKKSIYVF SKRKKRFLWF LHNFYVYECE YIFLFLRKQS SYLRSISSGV FLERTHFYGK IEYLRVMSCN SFQRILWFLK DIFIHYVRYQ GKAILVSKGT LILMNKWKFH FVNFWQFYFH FWFQPYRIRI KQLPNYSFSF MGYFSSVLKN PLVVRNQMLE NSFLINTLTN KLDTIAPVIF LIRLSKAQF CTVLGHPISK SIWTNLSDSD ILDRFCRICK NLCRYHSGSS KKQVLYRIKY ILRLSCARTL ARKHKSTVRT FMRRLGSAFL EEFFEEEEQS LSLVFLQKIP FLLHGLNRER IWYLDIIRMN DLVDHS
Source	Mammalian cell
Target Names	matK
Protein Names	Recommended name: Maturase K Alternative name(s): Intron maturase
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.