



Recombinant Arabidopsis thaliana Putative F-box/FBD/LRR-repeat protein At5g25850 (At5g25850)

Product Code	CSB-BP661659DOA
Abbreviation	At5g25850
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3E944
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MMGRKNSKQV RYSNRLILIN RLSQLSDPLI CQILSHLPIK EVVTTSVLST RWKNIWLSVP SLELIYSIFP NFNTFRSFGD RFFDSTRVSC IHNLKLYFDE HDACYLTSWI NAFVTRNIRR LYVRRVRGNY FHELPLSLYV CETLVSLKLF HLTLVDPEFV SLPCLKIMHL NYVWFPNDAT FERLVSSCPV LEDLKIDVLW NDGRVYRVKS RSLKRLRLLR SSTMVFDSPV GVVIDAPLIC YLRINDRVSE TYILNDLESN AKLDIHLDFG SEDFGEPSVS SRRSRIRSFL PAISKVMGLK ISENTFKVIH HYSNLAQLPQ FDNTSRLHVT LRVSELKWLP NFLECCPNLK SLIVAFNGDF EKMGSSEMNQ ISFLRMKRI VDMTIFMYTS KVLLHSDNIN QLSFSSVPKC LLSSLQFVEL NAQILRFDGE ILNLAKYFLE NSSILQKLT HPNKYGSTYA NMLHKFRRRS RTCQFILTDL E
Source	Baculovirus
Target Names	At5g25850
Protein Names	Recommended name: Putative F-box/FBD/LRR-repeat protein At5g25850
Expression Region	1-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.