



Recombinant Bubalus carabanensis Interleukin-1 alpha (IL1A)

Product Code	CSB-MP661675BAAW
Abbreviation	IL1A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3HWU1
Product Type	Recombinant Protein
Immunogen Species	Bubalus carabanensis (Swamp type water buffalo) (Bubalus bubalis carabanensis)
Purity	>85% (SDS-PAGE)
Sequence	SAHYSFQS NVKYNLLRVI HQECILNDAL NQSLIRDLSG PYLTAATLNN LEEAVKFDMMV AYWSEEDSQL PVTLRISKTKQ LRVSAQNEDE PVLLKEMPET PKIIKDETNL LFFWEKHGSM DYFKSVAHPK LFIATKQEKL VHMASGPPSV TDFQILEK
Source	Mammalian cell
Target Names	IL1A
Protein Names	Recommended name: Interleukin-1 alpha Short name= IL-1 alpha
Expression Region	113-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.