



Recombinant Chlamydia trachomatis serovar A 10 kDa chaperonin (groS)

Product Code	CSB-EP662043CAAZ-B
Abbreviation	groS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3KMQ8
Product Type	Recombinant Protein
Immunogen Species	Chlamydia trachomatis serovar A (strain ATCC VR-571B / DSM 19440 / HAR-13)
Purity	>85% (SDS-PAGE)
Sequence	MSDQATTLKI KPLGDRILVK REEEASTARG GIILPDTAKK KQDRAEVLAL GTGKKDDKGQ QLPFEVQVGD IVLIDKYSQG ELTVEGEEYV IVQMSEVIAV LQ
Source	E.coli
Target Names	groS
Protein Names	Recommended name: 10 kDa chaperonin Alternative name(s): 11.2 kDa stress response protein GroES protein Heat shock protein 10 Short name= HSP10 Protein Cpn10
Expression Region	1-102
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.