




Recombinant Human Innate immunity activator protein (INAVA)

Product Code	CSB-YP662049HU
Abbreviation	Recombinant Human C1orf106 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3KP66
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MESKDEVSDTDSGIILQSGPDSPVSPMKELTHAVHKQQRLEARLEACLEELR RLCLREAELTGTLPAEYPLKPGEKAPKVRRRRIGAAYKLDDWALHREDPLSSLE RQLALQLQITEAARRLCLEENLSRQARRQRKHSMLQEEKKLQELQRCLVERR RNSEPPPAALPLGRELSASDDSSLSDGLLLEEEESQVPKPPPEPAPPSPRPL PPQTLEGLQPTGPEAGSPERAPVQNSPWKETSLDHPYEKPRKSSEPWSESS SPATTPQDGPSASSLWLEPASYHVVPPIRGVPGQWQGRTSAPATPEIQGRRG QSQSLRVDSFRAGPEGRGRSAFPRRRPTHYTVTVPDSCFPATKPLPHAACH SCSEDSGSDVSSISHPTSPGSSSPDISFLQPLSPPKTHRRGAWVPAGSRELV AHHPKLLPPGYFPAGRYVVVAESPLPPGEWELRRAAPGPAYEEEGTPLYRQ RLVPSRSRIVRTPSLKDSPAGRGLSKAAVSEELKWWHERARLRSTRPHSLDR QGAFRVRSPLPLGREGFGRALGPRAQVPTVCVLRRLSPDGAPVQVFVPEKGEIIS QV
Research Area	Cell Biology
Source	Yeast
Target Names	INAVA
Protein Names	Recommended name: Uncharacterized protein C1orf106
Expression Region	1-578aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	65.5kDa
Protein Length	Full Length of Isoform 2
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.