



Recombinant Axin-like protein 1 (axl-1)

Product Code	CSB-EP662072CXY
Abbreviation	axl-1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q3LRZ3
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTMRSKFSID RVLHDTAF AK WCQKEESRNA DSITLYTSIL DFESKLKTGS PSLSLLKLAR HIHRKYVSLN TGTCHVIGDS IRTKISTRVH EVLDGKPPYI DLFDPLKQPL FQHLRSMHTE FSTTTADVNT TWEDASSTSS SNKGATIWFN EDAIDKSSRH EIGQSTVTHE SEDDRFSFFN AVCTRLNSLQ ETKNSSETEE HAESPRKEKS STPYGTDGFA PPRSTQTNT LKVSNLPKRF ESLYKKKRQQ VATSDSSGFG SNASDFWSFE RYGKSNQGT L ERPNRLFTGT NNGFSTLQPK RRGSEAPKMT VELRYENDVP MVAKISANHA QSVTLRYFRH LFGLHYTDNC RFFFKSTCED GSAQYQWTL L FHDDDILPVF QNRITAICRM CPPPEDHDLI
Source	E.coli
Target Names	axl-1
Protein Names	Recommended name: Axin-like protein 1
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>