



# Recombinant Bovine Na (+)/H (+) exchange regulatory cofactor NHE-RF1 (SLC9A3R1)

<b>Product Code</b>	CSB-MP662264BO
<b>Abbreviation</b>	SLC9A3R1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3SZK8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SADAGAGAP LPRLCCLEKG PNGYGFHLHG EKGKVGQYIR LVEPGSPAЕК SGLLAGDRLV EVNGENVEKE THQQVVNRIR AALNSVRLLV VDPETDERLQ KLG VQVREEM LRAQEGPGQA EPPAAAAEER GAGGENEPPA AAPEPREAEQ SPQERREL RP RLCAMKKGPN GYGFNLHSDK SKPGQFIRAV DPDSPAEASG LRAQDRIVEV NGVCVEGKPH GEVVSАIKAG GDEAKLLVVD RETDEFFKKC KVIPSQEHLQ GPLPEPITNG EIEKENSPEA LAETASESPM PPLARTASSD TSEELNSQDS PPKQDSTAPS STSSSSSDPV LDFSISLAVA KERAHQKRVS KRAPQMDWSK KNELFSNL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SLC9A3R1
<b>Protein Names</b>	Recommended name: Na(+)/H(+) exchange regulatory cofactor NHE-RF1 Short name= NHERF-1 Alternative name(s): Ezrin-radixin-moesin-binding phosphoprotein 50 Short name= EBP50 Regulatory cofactor of Na(+)/H(+) exchanger Sodium-hy
<b>Expression Region</b>	2-368
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.