



Recombinant Bovine Pro-cathepsin H (CTSH)

Product Code	CSB-EP662282BO-B
Abbreviation	CTSH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3T0I2
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	YPPSM DWRKKGNFVT PVKNQGSCGS CWFSTT GAL ESAVAIATGK LPFLAEQQLV DCAQNFNNHG CQGGLPSQAF EYIRYNKGIM GEDTYPYRGG DGDCKYQPSK AIAFVKDVAN ITLNDEEAMV EAVALHNPVS FAFEVTADFM MYRKGIYSST SCHKTPDKVN HAVLAVGYGE EKGIPYWIVK NSWGPNWGMK GYFLIERGKN MCGLAACASF PIPLV
Source	E.coli
Target Names	CTSH
Protein Names	Recommended name: Pro-cathepsin H Cleaved into the following 4 chains: 1. Cathepsin H mini chain 2. Cathepsin H EC= 3. 3.4.22.16 4. Cathepsin H heavy chain 5. Cathepsin H light chain
Expression Region	116-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a lysosomal cysteine proteinase important in the overall degradation of lysosomal proteins. It is composed of a dimer of disulfide-linked heavy and light chains, both produced from a single protein precursor. The encoded protein, which belongs to the peptidase C1 protein family, can act both as an aminopeptidase and as an endopeptidase. Increased expression of this gene has been correlated with malignant progression of prostate tumors. Two transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.