



Recombinant Mouse F-box/WD repeat-containing protein 1A (Btrc)

Product Code	CSB-MP662374MO
Abbreviation	Btrc
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3ULA2
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MDPAEAVLQE KALFKMCSMP RSLWLGCSL ADSMPSLRCL YNPGTGALTA FQNSSEREDC NNGEPPRKII PEKNSLRQTY NSCARLCINQ ETVCLTSTAM KTENCVAKAK LANGTSSMIV PKQRKLSASY EKEKELCVKY FEQWSESDQV EFVEHLISQM CHYQHGHINS YLKPMLQRDF ITALPARGLD HIAENILSYL DAKSLCAAEL VCKEWYRVTS DGMLWKKLIE RMVRTDSLWR GLAERRGWGQ YLFKNKPPDE NAPPNSFYRA LYPKIIQDIE TIESNWRGCR HSLQRIHCRS ETSKGVYCLQ YDDQKIVSGL RDNTIKIWDK STLECKRILT GHTGSVLCLQ YDERVIITGS SDSTVRVWDV NAGEMLNTLI HHCEAVLHLR FNNGMMVTCS KDRSIAVWDM ASPTDITLRR VLVGHRAAVN VVDFDDKYIV SASGDRTIKV WNTSTCEFVR TLNGHKRGIA CLQYRDRLVV SGSSDNTIRL WDIECGACLR VLEGHEELVR CIRFDNKRIV SGAYDGKIKV WDLMAALDPR APAGTLCLRT LVEHSGRVFR LQFDEFQIVS SSHDDTILIW DFLNDPAAHA EPPRSPSRTY TYISR
Source	Mammalian cell
Target Names	Btrc
Protein Names	Recommended name: F-box/WD repeat-containing protein 1A Alternative name(s): Beta-TrCP protein E3RS-IkappaB Beta-transducin repeat-containing protein Short name= Beta-TrCP E3RS-IkappaB Short name= mE3RS-IkappaB F-box and W
Expression Region	1-605
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.