



# Recombinant Arabidopsis thaliana Probable protein phosphatase 2C 49 (At3g62260)

<b>Product Code</b>	CSB-BP662983DOA
<b>Abbreviation</b>	At3g62260
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3EAF9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVAEAEVVFQ QSLPAVLEIE LFDGVSSAVK SPVSSPKLGF TQSTASVSGS LTTSPVQADI FPEGDCPSV LDYIPTIRSG SFADIGPKRN MEDEHIRIDD LSSQVGSLEF LPKPSAFYAV FDGHGGPEAA AYPVRENAIRF FFEDEQFPQT SEVSSVYVEE VETSLRNAFL QADLALAEDC SISDSCGTTA LTALICGRLL MVANAGDCRA VLCKRGRAID MSEDHKPINL LERRRVEESG GFITNDGYLN EVLAVTRALG DWDLKLPHGS QSPLISEPEI KQITLTEDDE FLVIGCDGIW DVLTSQEAVS IVRRGLNRHN DPTRCARELV MEALGRNSFD NLTAVVVCFM TMDRGDKPVV PLEKRRCFSL SPEAFCSLRN LLDG
<b>Source</b>	Baculovirus
<b>Target Names</b>	At3g62260
<b>Protein Names</b>	Recommended name: Probable protein phosphatase 2C 49 Short name= AtPP2C49 EC= 3.1.3.16
<b>Expression Region</b>	1-384
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.