



Recombinant Arabidopsis thaliana Putative F-box protein At1g53360 (At1g53360)

Product Code	CSB-MP662989DOA
Abbreviation	At1g53360
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3ECR3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNQSLRRGY WFWFKYEPNL FWYYSIPTKH ANIILVYLLI EMFSRVPAKS IARFRCVSES IHRRPDVTEQ SLTKSMSRPR LLFALEVNKD LLFFSSPQPQ NPYDNSSLVA TPYKRFPKYL PTRICTTLGG LVFLQKWLRK KTRVICNPVS GDYITLPKVK ATGVGESYFG FDPITKQFKV LCMTWSRYGT PNTHQVLTLE TGKRLWRTIQ DPILPHYRSF DRICINGVLY YGADFEEESQS SKIVCFDFRF EKFSFINIAD EGMFRGSYKW TLFNYKGKLG AHQYSRNGEL VLWVLEDAEN HKWSKSICIL PPIVHNSRIV GVTGTGEIVF SPYACYMPSP FYIFFYNIQT KTCTRVHVKG FEEFKHNFTL LHTFLDFVED TKFM
Source	Mammalian cell
Target Names	At1g53360
Protein Names	Recommended name: Putative F-box protein At1g53360
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.