



Recombinant Natronomonas pharaonis Signal recognition particle 19 kDa protein (srp19)

Product Code	CSB-MP663065NAAJ
Abbreviation	srp19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3IPF4
Product Type	Recombinant Protein
Immunogen Species	Natronomonas pharaonis (strain ATCC 35678 / DSM 2160 / CIP 103997 / NBRC 14720 / NCIMB 2260 / Gabara) (Halobacterium pharaonis)
Purity	>85% (SDS-PAGE)
Sequence	MVENVIWPAY LDANCTRNEG RRVAEDLAVE EPTVDEIASA VQQVGYDAVI ERSKTYPREY EPRGRVVVKG ADDASKSDLL EAVAAYVTAI RE
Source	Mammalian cell
Target Names	srp19
Protein Names	Recommended name: Signal recognition particle 19 kDa protein Short name=SRP19
Expression Region	1-92
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.