



Recombinant Bovine Exosome complex component RRP40 (EXOSC3)

Product Code	CSB-MP663601BO
Abbreviation	EXOSC3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3T0E1
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	AEAAGVPAE SLAGCRARAA RTVLDQVVLP GEELLLPDQE DGDGPGGAGE RPLRLNAAAR SRGRVVC GPG LRRCGDRLLV TKCGRLRHKE PGSGSGGGVY WVDSQQKRYV PVKGDHVI GVTAKSGDIFK VDVGGSEPAS LSYLA FEGAT KRNRPNVQVG DLIYGQFVVA NKDMEPEMVC IDSCGRANGM GVIGQDGLLF KVTGLIRKL LAPDCEILQE VGKLYPLEIV FGMNGRIWVK AKTIQQTLIL ANILEACEHM TADQRKQIFS RLAES
Source	Mammalian cell
Target Names	EXOSC3
Protein Names	Recommended name: Exosome complex component RRP40 Alternative name(s): Exosome component 3 Ribosomal RNA-processing protein 40
Expression Region	2-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.