



# Recombinant Mouse CREB-regulated transcription coactivator 2 (Crtc2)

<b>Product Code</b>	CSB-EP663643MO
<b>Abbreviation</b>	Crtc2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3U182
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ATSGANGPG SATASASNPR KFSEKIALQK QRQAEETA AF EEVMMDIGST RLQAQKLRLA YTRSSHYGGS LPNVNQIGCG LAEFQSPLHS PLDSSRSTRH HGLVERVQRD ARRMVSPRR YPRHIDSSPY SPAYLSPPE SGWRRMMPWG NLPAEKGQLF RLPSALNRTS SDSALHTSVM NPNPQDTYPG PTPPSVLPSR RGGGFLDGEM DAKVPAIEEN VVDDKHLLKP WDAKKLSSSS SRPRSCEVPG INIFSPDQP ANVPVLPPAM NTGGSLPDLT NLHFPPPLPT PLDPEETVYP SLGGNSTTN LHTMTHLGI SGGLGLGPSY DVPGLHSPLS HPSLQSSLSN PNLQASLSSP QPQLQGSHSH PSLPASSLAH HALPTTSLGH PLSAPALSS SSSSSSTSSP VLSAPPYPAS TPGASPRHRR VPLSPLSLPA GPADARRSQQ QLPKQFSPTM SPTLSSITQG VPLDTSKLPT DQRLPPYPYS PPSLVIPTHP PTPKSLQQLP SQACLVQPSG GQPPGRQPHY GALYPPGSSG HGQQPYHRPI NDFSLGNLEQ FNMESPSTSL VLDPPAFSEG PGFLGSEGSM SGPQDPHVLN HQNLTHCSRH GSGPNIILTG DSSPGFSKEI AAALAGVPGF EVSASGLELG LGLEDEL RME PLGLEGLTML SDPCALLPDP AVEDSFRSDR LQ
<b>Source</b>	E.coli
<b>Target Names</b>	Crtc2
<b>Protein Names</b>	Recommended name: CREB-regulated transcription coactivator 2 Alternative name(s): Transducer of regulated cAMP response element-binding protein 2 Short name= TORC-2 Short name= Transducer of CREB protein 2
<b>Expression Region</b>	2-692
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.