



Recombinant Arabidopsis thaliana Putative F-box/LRR-repeat protein At5g25860 (At5g25860)

Product Code	CSB-BP664252DOA
Abbreviation	At5g25860
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3E7J7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDTRKVKSSS RDAVNCLPDE ILAKILSYLP TKRAVSTSLI SKRWRNLFAL MIQLFESQHH LYLDDSDLVY PEEGKGEMKD VQESFGDFVD KTLTDCNTIK KLSILCPLKC CAHKDIDQWL HHAMERGVVD LDMRFKMGFT PTQNEWPCNV FTIKTLVKLT FRIEVGQNKM PYCPKGVLPV LKSLFLHAIW FACERLCHSM LPGCPILEEL FLHDLRGYRY NDSPNFSISH KTLKRLTVHF NNSFELGRPM QFDTPSLLYL EYSGFAPCRY SRTSYLNSLV EVKLDVHIRH ECNTHMHLSP IIDWIGNVKT LSLSPASVKM MYSRCVELHV FSNLVKLYFE SNEKEGWVFL PRLLNKSPKL ETLVLKGLHC ASDHGCVDNR NEVKVLEIYG FSGCGREVRQ VKCLLREMQF LQVMKVEIDA HDDNKCLRVI NHLLDLPKRP SKFHIQFL
Source	Baculovirus
Target Names	At5g25860
Protein Names	Recommended name: Putative F-box/LRR-repeat protein At5g25860
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.