



# Recombinant Arabidopsis thaliana Putative F-box/LRR-repeat protein At5g25860 (At5g25860)

<b>Product Code</b>	CSB-MP664252DOA
<b>Abbreviation</b>	At5g25860
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3E7J7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDTRKVKSSS RDAVNCLPDE ILAKILSYLP TKRAVSTSLI SKRWRNLFAL MIQLFESQHH LYLDDSDLVY PEEGKGEMKD VQESFGDFVD KTLTDCNTIK KLSILCPLKC CAHKDIDQWL HHAMERGVVD LDMRFKMGFT PTQNEWPCNV FTIKTLVKLT FRIEVGQNKM PYCPKGVLPV LKSLFLHAIW FACERLCHSM LPGCPILEEL FLHDLRGYRY NDSPNFSISH KTLKRLTVHF NNSFELGRPM QFDTPSLLYL EYSGFAPCRY SRTSYLNSLV EVKLDVHIRH ECNTHMHLSP IIDWIGNVKT LSLSPASVKM MYSRCVELHV FSNLVKLYFE SNEKEGWVFL PRLLNKSPKL ETLVLKGLHC ASDHGVCVDR NEVKVLEIYG FSGCGREVRQ VKCLLREMQF LQVMKVEIDA HDDNKCLRVI NHLLDLPKRP SKFHIQFL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At5g25860
<b>Protein Names</b>	Recommended name: Putative F-box/LRR-repeat protein At5g25860
<b>Expression Region</b>	1-448
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.